



www.IetLtd.com Proudly serving laboratories worldwide since 1979

CALL +1.847.913.0777 for Refurbished & Certified Lab Equipment

Perkin Elmer Series 200 LC Pumps

Isocratic, Binary, Quaternary and Micro Flow
Versions Built to Withstand the Test of Time

Top performance characteristics are what you have come to expect from PerkinElmer™ LC pumps. We've taken this exceptional patented design and put it all into a common, upgradable pump platform. The PerkinElmer Series 200 Pumps - a fifth-generation pump with superior flow rate precision, compositional accuracy and exceptional retention time repeatability for all of your separations - isocratic or gradient. While the new Series 200 Micro Pumps provide precise high pressure gradients and pulseless flow for microbore, narrowbore, LC/MS and analytical requirements.

The PerkinElmer Series 200 Pump provides stand-alone capability or full external control through TotalChrom™ Chromatography Workstation. Start out with the isocratic version now and upgrade to full binary or quaternary gradient capability - all within the same pump - we'll ship you a cost-effective upgrade kit for easy installation right in your lab.

Proven Rugged & Reliable Design

Upgradable isocratic, binary, and quaternary versions as well as biocompatible options for maximum flexibility.

- Micro Pump format for ultra-low flow rate applications.
- Small, compact design [12 inches (30.5 cm) wide] saves valuable bench space.
- Large back-lit LCD screen allows clear viewing of pump method and status parameters at all times.
- Solvent resistant front panel with rounded ergonomic design for easier interaction with keypad.
- TotalChrom ready — just connect to the link interface and go! No upgrades or additional cables required!
- Tactile keyboard with unique “key click” so that you know your entry was made.
- Snap-in contact closure connections and finger-tight plumbing fittings make it simple to connect or disconnect your system.

- Visual access front-door panel allows quick plumbing inspection and easy access to prime purge valve.
- Pump module can be accessed or completely removed to lab bench for maintenance without having to “unstack” modules via unique slide-out design.

Comprehensive Method and Editing Power

- Directory of up to 20 methods, 20 steps each with built-in method linking increases sample throughput.
- Edit the active method in real-time with easy access function keys.
- Automatic shutdown of the pump method at your final gradient conditions.
- Graphic display of solvent profile on-screen shows real-time modifications to gradient method.
- Hold gradients or advance to the next step with a simple keystroke.

Complete Diagnostics and Method Documentation Simplifies GLP Compliance

- Full method documentation including modifications with method name, number, date and time stamp.
- Maintenance Log coded by the serial number of the pump, including pump and seal cycles.
- Extensive service diagnostics.
- Instrument Log records history of operation for quick diagnosis and servicing to get you up and running in no time.

PerkinElmer’s Pump Design Offers the Best Flow Rate Precision. Here’s Why...

A Automated Solvent Compensation Maintains a continuous and constant flow with each pump stroke using patented real-time blending compensation, which accounts for and meters different volumes according to the solvent shrinkage that occurs with mixing. This same principle corrects for solvent compressibility changes.

B Minimal Pulsations Minimizes pressure pulsations due to high speed retraction of delivery piston with every stroke (65 msec.)

C Low Mixing Volume* Provides high-efficiency solvent mixing resulting from extremely high turbulence inside the pump when the pump transfers the next increment of solvent from the low-pressure piston to the high-pressure piston at a high flow rate (100 mL/min).

D No Air Bubbles Eliminates cavitation or bubbles as a result of the unique, patented design that provides positive transfer of the next increment of solvent under pressure rather than vacuum refilling, as is typically done with other pumps. In the schematic, the delivery chamber is about to be positively refilled from the metering chamber.

E Superior Low Flow-Rate Performance Operates efficiently even at lower flow rates because of the double check valve design that prevents solvent backflow on each delivery stroke.

Time-Proven Design Keeps Your Lab Going

PerkinElmer pumps are built to keep running. We ensure that the pump going into your laboratory is fully tested and calibrated to ensure reliable, trouble-free operation all day, every day. Our quality assurance laboratory reports that the typical Mean Time Between Failure (MTBF) is a minimum of 4400 hours - that's two years of continuous operation with no failures of any kind (8 hours a day, 5 days a week). PerkinElmer pumps are still the lowest maintenance pumps on the market with the most rugged design.

Time-Proven Design Provides Outstanding Flow Rate Reproducibility, Superior Compositional Accuracy and Unsurpassed Gradient Performance

- *Unsurpassed Compositional Accuracy.* The proportioning of the Series 200 Pump has been designed so that the minor component is always “sandwiched” between the major components. In this way, compositional accuracy is preserved even when one mobile phase component is mixed at the very low percentage.

- *No External Mixing Volume Is Required**. All mixing is done inside the pump. This is accomplished after the metering stroke when proportioned solvent is transferred to the high pressure side at 100mL/min. The resulting turbulence causes thorough mixing and eliminates the need for an external dynamic mixer.

Series 200 Pump Specifications

	Series 200 Micro Pumps	Series 200 Analytical Pumps
Flow Rate Range	1µL/min - 3mL/min	0.01 to 10.0 mL/min
Flow Rate Increments	1µL/min from 0 to 0.999 mL/min; 0.01 to 3 mL/min	0.01 from 0 to 0.99 mL/min; 0.1 from 1.0 to 10 mL/min
Flow Precision	0.3% RSD (typical 0.1%) at 300µL/min water @ 1000 psi	0.3% RSD (typical 0.1%) at 1 mL/min water @ 1000 psi
Flow Accuracy	+/- 1% of setting at 300µL/min and 1000 psi with water	+/- 1% of setting at 1 mL/min and 1000 psi with water
Retention Time Reproducibility	<0.3% RSD (typically <0.1%)	<0.3% RSD (typically <0.1%)
Pressure Range	0 to 6200 psi for flow-rate range	0 to 6200 psi for flow-rate range
Compositional Range	0 to 100%, settable to 1%	0 to 100%, solvent A to D, settable to 0.1%
Compositional Accuracy	Typically 0.5% from 10 to 90% up to 2mL/min	Typically 0.5% from 3 to 97% up to 5mL/min
Compositional Precision Gradient Forms	Typically < 0.2% variation Linear, convex and concave w/ high pressure blender	Typically < 0.2% variation Linear, exponential 1 to 9.9 (positive and negative) or step (89 profile selections)
Linearity	Within 1% from 10 to 90%	
Program Steps	Step 0 (Equil) plus up to 20 solvent program steps. Battery-backed storage of up to 20 methods plus 1 edit method and 1 default method	
Step Time	0 to 999 minutes	
Step Time	Increments 0.1 min from 0 to 9.9 min; 1.0 min from 10 to 999 min	

Timed Events
Events / Method
Event Increments

Two built-in independent momentary (2 second) contact closures
Max of 9
0.1 min from 0 to 9.9 min; 1.0 min from 10 to 999 min

Part Number Description

N291-0100 Series 200 LC Pump Isocratic Version (50/60 Hz)
N291-0101 Series 200 LC Pump Binary Version (50/60 Hz)
N291-0102 Series 200 LC Pump Quaternary Version (50/60 Hz)
N291-0202 Series 200 BIO LC Pump Quaternary Version (50/60 Hz)
N291-0500 Series 200 Micro Pump Isocratic Version (50/60 Hz)
N291-0503 Series 200 Micro Pump Binary Package (50/60 Hz) (includes a Series 200 Dynamic Mixer)
N291-0342 Series 200 Upgrade Kit – Isocratic to Binary
N291-0343 Series 200 Upgrade Kit – Isocratic to Quaternary
N291-0344 Series 200 Upgrade Kit – Binary to Quaternary



Proudly serving laboratories worldwide since 1979

CALL +1.847.913.0777 for Refurbished & Certified Lab Equipment