



[www.letLtd.com](http://www.letLtd.com) Proudly serving laboratories worldwide since 1979

**CALL 001.847.913.0777 for Refurbished & Certified Lab Equipment**

## ***Typhoon 9410 Gel and Blot Imager Data Sheet***

Typhoon Variable Mode Imager

Multicolor fluorescence, filmless autoradiography, and chemiluminescence

Amersham Biosciences has been manufacturing highly sensitive imaging instrumentation for over a decade. Typhoon™ Variable Mode Imager unites proven storage phosphor autoradiography technology with four-color fluorescent labelling techniques for maximum data quality in a single, high-throughput system (Fig 1).

For DNA, RNA, and protein samples, users can choose from:

- storage phosphor autoradiography
- direct green (532 nm) excited fluorescence
- direct red (633 nm) excited fluorescence
- direct blue (457 nm) excited fluorescence (Typhoon 9400/9410 only)
- direct blue (488 nm) excited fluorescence (Typhoon 9400/9410 only)
- chemiluminescence

When one of the six scanning modes is selected, the appropriate optical components are automatically activated. Typhoon scans mounted and unmounted storage phosphor screens, plus gels and blots up to 35 × 43 cm as well as microarray slides for Typhoon 9210 and 9410 models. Typhoon Variable Mode Imager exhibits outstanding linearity and quantitative accuracy, and includes ImageQuant™ Image Analysis Software for Windows™ 2000.

Typhoon is fully optimized as part of the Ettan™ DIGE system, and seamless integration with DeCyder™ Differential Analysis Software is ensured for all models. Special tray templates that are part of the latest scanner control software and new optional Gel Alignment Guides add higher throughput and ease of handling 2-D gels run on Ettan DALT and SE600.

| <b>Typhoon Model</b> | <b>9200</b> | <b>9210</b> | <b>9400</b> | <b>9410</b> |
|----------------------|-------------|-------------|-------------|-------------|
| Phosphorimaging      | Yes         | Yes         | Yes         | Yes         |
| 633 nm-Fluorescence  | Yes         | Yes         | Yes         | Yes         |
| 532 nm-Fluorescence  | Yes         | Yes         | Yes         | Yes         |
| 488 nm-Fluorescence  | —           | —           | Yes         | Yes         |
| 457 nm-Fluorescence  | —           | —           | Yes         | Yes         |
| Chemiluminescence    | Yes         | Yes         | Yes         | Yes         |
| Microarray           | —           | Yes         | —           | Yes         |

**Table 1.** Each Typhoon model offers different features, and for maximum flexibility, each model has an upgrade path.

### **Components**

- Typhoon scanner with TCP/IP, scan control software for Windows 2000
- ImageQuant Image Analysis Software for Windows 2000
- User's guide
- Microarray slide holder (Typhoon 9210/9410 only)
- Ettan DIGE components sold separately

## **Detection threshold**

### **Storage phosphor**

Storage phosphor screens retain energy from beta particles, X-rays, and gamma rays. The lower limit of detection for a 1-h exposure is less than 2 dpm/mm<sup>2</sup> for <sup>14</sup>C (200 µm only). The lower limit of detection for <sup>32</sup>P is typically 5–10 times lower than the limit for <sup>14</sup>C.

- 488 nm-excited fluorescence (Typhoon 9400/9410 only) 100 amol fluorescein end-labelled DNA oligonucleotide in 12% polyacrylamide gel sandwich, 0.4 mm thick.
- 532 nm-excited fluorescence 200 amol HEX<sup>TM</sup>, TAMRA<sup>TM</sup>, ROX<sup>TM</sup>, and 400 amol fluorescein end-labelled DNA oligo-nucleotide in 12% polyacrylamide gel sandwich, 0.4 mm thick.
- 633 nm-excited fluorescence 200 amol Cy<sup>TM</sup>5 end-labelled DNA oligonucleotide in 12% polyacrylamide gel sandwich, 0.4 mm thick.

## **Light measurement**

### **Storage phosphor**

Light is emitted from the storage phosphor screen in proportion to the amount of radioactivity in the sample upon laser-induced stimulation. Emitted light is collected and converted to an electrical signal by a photomultiplier. The electrical signal is digitized to permit image display and analysis.

### **Fluorescence**

Upon excitation, light is emitted from a fluorescently-labelled sample in proportion to the amount of labelled compound in the sample. Emitted light is collected and converted to an electrical signal by a photomultiplier. The electrical signal is digitized for image display and analysis.

### **Chemiluminescence**

Emitted light from a chemiluminescent reaction is collected and converted to an electrical signal by a photomultiplier. The electrical signal is digitized for image display and analysis.

## Data Storage

Data are stored in a square root encoded 16-bit TIFF to provide the digital resolution required to characterize subtle signal intensity differences over the wide dynamic range of the instrument.

## Specifications

### Spatial Resolution

|                                  | <b>Typhoon 9200</b> | <b>Typhoon 9210</b> | <b>Typhoon 9400</b> | <b>Typhoon 9410</b> |
|----------------------------------|---------------------|---------------------|---------------------|---------------------|
| Storage phosphor autoradiography | 2 line pairs/mm     | 2 line pairs/mm     | 2 line pairs/mm     | 2 line pairs/mm     |
| Green-excited fluorescence       | 8 line pairs/mm     | 10 line pairs/mm    | 8 line pairs/mm     | 10 line pairs/mm    |
| Red-excited fluorescence         | 8 line pairs/mm     | 10 line pairs/mm    | 8 line pairs/mm     | 10 line pairs/mm    |
| Blue-excited fluorescence        | N/a                 | N/a                 | 8 line pairs/mm     | 10 line pairs/mm    |

### Emission filters

| <b>Filter type</b>           | <b>Wavelength range (nm)</b> | <b>Fluorochromes</b>                        |
|------------------------------|------------------------------|---------------------------------------------|
| 520 nm-bandpass (520 BP 40)* | 500–540                      | Cy2, ECL Plus™, Vistra Green™ fluorescein** |
| 555 nm-bandpass (555 BP 20)  | 545–565                      | R6G, HEX                                    |
| 580 nm-bandpass (580 BP 30)  | 565–595                      | Cy3, TAMRA                                  |
| 610 nm-bandpass (610 BP 30)  | 595–625                      | ROX, ethidium bromide, SYPRO™ Ruby          |
| 670 nm-bandpass (670 BP 30)  | 655–685                      | Cy5                                         |
| 526 nm-short-pass (526 SP)   | < = 526                      | Fluorescein\$, Vistra Green                 |
| 560 nm-long-pass (560 LP)    | <= 560                       | TRITC                                       |

\* For Typhoon 9400/9410 models only.

\*\* With blue excitation

\$ With green excitation

### Exposure time

Typically, storage phosphor screen exposure takes 10% of the time for an equivalent exposure to conventional film.

### Single or dual-channel scanning time (min)

| pixel size (µm) | 1000* | 500 | 200 | 100 | 50 | 25  |
|-----------------|-------|-----|-----|-----|----|-----|
| 20 × 25 cm      | 1     | 2   | 5   | 9   | 19 | 92  |
| 35 × 43 cm      | 2     | 4   | 10  | 21  | 40 | 167 |

*For Typhoon 9210/9410 models only*

|                      |    |
|----------------------|----|
| pixel size (µm)      | 10 |
| 2 × 8 cm (1 slide)   | 8  |
| 2 × 18 cm (2 slides) | 16 |

### Green light source

- Type: 20 mW solid state, doubled frequency SYAG laser
- Estimated average lifetime: ~10 000 h
- Wavelength: 532 nm

**Blue light source (Typhoon 9400/9410 only)**

- Type: 30 mW, Argon ion laser
- Estimated average lifetime: ~5 000 h
- Wavelength: 488 nm (20 mW), 457 nm (4 mW)

**Four-channel linked scanning time (min)**

|                 |   |       |     |     |     |     |    |
|-----------------|---|-------|-----|-----|-----|-----|----|
| pixel size (µm) |   | 1000* | 500 | 200 | 100 | 50  | 25 |
| 20 × 25 cm      | 3 | 5     | 9   | 19  | 37  | 186 |    |
| 35 × 43 cm      | 6 | 10    | 21  | 40  | 80  | 335 |    |

*For Typhoon 9210/9410 models only*

|                      |    |
|----------------------|----|
| pixel size (µm)      | 10 |
| 2 × 8 cm (1 slide)   | 17 |
| 2 × 18 cm (2 slides) | 32 |

*\*Not recommended for quantitative analysis.*

**Uniformity**

± 5% over entire scan area

**Pixel accuracy**

± 0.15%

**Data format**

16-bit (65 536 levels), TIFF (.GEL file extension)

**Linearity**

Less than 7.5% relative standard deviation for entire dynamic range

**Linear dynamic range**

Five orders of magnitude (100000:1)

**External interface**

10 Base-T Ethernet using the TCP/IP protocol

**Software**

- Scan control software for Windows 2000
- ImageQuant image analysis software for Windows 2000

**Red light source**

- Type: 10 mW Helium-Neon laser
- Estimated average lifetime: ~10 000 h
- Wavelength: 632.8 nm



[www.ietLtd.com](http://www.ietLtd.com)

Proudly serving laboratories worldwide since 1979

**CALL 001.847.913.0777 for Refurbished & Certified Lab Equipment**