



[www.IetLtd.com](http://www.IetLtd.com) Proudly serving laboratories worldwide since 1979

CALL +847.913.0777 for Refurbished & Certified Lab Equipment

## ***G1330B Autosampler Thermostat***

### **Performance Specifications Agilent 1200 autosampler thermostat**

Type	Specification
Temperature range:	setable from 4°C to 40°C in 1° increments
Temperature accuracy at ambient temperatures < 25°C and humidity < 50%	-1°C to +4°C at a set point of 4°C
Temperature accuracy at ambient temperatures > 25°C and/or humidity > 50%	-1°C to +5°C at a set point of 4°C

### **Physical Specifications - Thermostatted Autosampler**

Type	Specification	Comments
Thermostat Weight	20.7 kg (45.6 lbs)	
Dimensions (Height x width x depth)	140 x 345 x 435 mm (5.5 x 13.5 x 17 inches)	
Line voltage	100 – 120 or 220 – 240 VAC, ± 10 %	Automatic selection
Line frequency	50 or 60 Hz, ± 5 %	
Power consumption	Autosampler: 300 VA ALS Thermostat: 260 VA	Maximum Maximum
Ambient operating temperature	4 – 40 °C (41 – 131 °F)	See User manual
Ambient non-operating temperature	-40 – 70 °C (-4 – 158 °F)	
Humidity	< 95 %, at 25 – 40 °C (77 – 104 °F)	Non-condensing
Operating Altitude	Up to 2000 m (6500 ft)	
Non-operating altitude	Up to 4600 m (14950 ft)	For storing the autosampler
Safety standards: IEC, CSA, UL, EN	Installation Category II, Pollution Degree 2	



[www.letLtd.com](http://www.letLtd.com) Proudly serving laboratories worldwide since 1979

CALL +847.913.0777 for Refurbished & Certified Lab Equipment