



**IET | International Equipment Trading Ltd.**

[www.IetLtd.com](http://www.IetLtd.com) | *Proudly serving laboratories worldwide since 1979*

CALL +847.913.0777 for Refurbished & Certified Lab Equipment

## Jasco CO-2060/2065 Column Oven Performance Specifications

The CO-2060/2065 ovens are capable of temperature control with a precision of  $\pm 0.1^{\circ}\text{C}$  with excellent reproducibility. The CO-2060 is for heating and cooling, while the CO-2065 is for heating only. Both models detect abnormal temperatures and solvent leaks and include a function for shutting off the heat.

A manual indicator, preheating coil, reaction coil, and other such devices can be installed on both the CO-2060 and the CO-2065. Furthermore, the CO-2060 can easily perform not only cooling (used for optical isomer separation columns), but also temperature increases and decreases by the use of a time program.

### Specifications:

Model name	CO-2060	CO-2065
Temperature control method	Heater, Peltier	Heater
Temperature control range	15 °C below ambient to 80°C	10°C above ambient to 80°C 20°C above ambient to
Accuracy	$\pm 0.1^{\circ}\text{C}$	
	(When temperature is set at 40°C.)	
Stirrer	None	
Safety check	Heater power shut off when unusually high temperature or flammable solvent leaks are detected.	
Column compartment	110(W) × 415(H) × 83(D) mm	
Dimensions	150(W) × 465(H) × 470(D) mm	

Weight	Approx. 19 kg	Approx. 15 kg
Power requirements	100~220V, AC, 50/60Hz	
Power consumption	Approx. 355VA	Approx. 230VA



WWW.IETI.TD.COM



**IET | International Equipment Trading Ltd.**

[www.IetLtd.com](http://www.IetLtd.com) | *Proudly serving laboratories worldwide since 1979*

CALL +847.913.0777 for Refurbished & Certified Lab Equipment

**IET - Used Lab Equipment - Refurbished Analytical Laboratory Instruments**