

IET | International Equipment Trading Ltd. <u>www.IetLtd.com</u> | *Proudly serving laboratories worldwide since 1979* CALL +847.913.0777 for Refurbished & Certified Lab Equipment

Mettler XP 204

Performance Specifications

Features and Benefits:

SmartGrid: unparalleled weighing performance, thanks to the revolutionary grid weighing pan which successfully minimizes the effects of turbulence in the weighing chamber. Results are more accurate and stabilization times dramatically shorter.

SmartSens: hands-free operation for personnel security The user can concentrate fully on the safe handling of samples to avoid spillage and minimize the risk of contamination.

SmartScreen: guides the user by color for safe operation. Outstanding ease of use with touch-screen technology.

LevelControl: warns if the balance is not correctly levelled to ensure the accuracy of your results.

Connectivity: pace-setting interfacing flexibility - including Ethernet, Bluetooth (wireless connection) and PS/2 - for efficient data capture and easy network integration.

Integrated anti-static kit (option): discharge static samples instantly without any risk of contamination.

Specifications - XP204 Analytical Balance	
Maximum Capacity	220 g
Readability	0.1 mg
Minimum weigh (typical acc. USP)	120 mg
Minimum weigh (typical, U=1%, sd=2)	8 mg
Taring range	0220 g
Repeatability	0.05 - 0.07 mg
Linearity	0.2 mg

Eccentric Load	0.25 mg
Adjustment with internal weights	proFACT, fully automatic time- and temperature- controlled adjustment
Adjustment with external weights	custom weights
Sensitivity Accuracy	3.0 x 10^-6 •Rnt
Sensitivity Stability	1 x 10^-6/a •Rnt
Sensitivity offset	3x10-6·Rnt
Sensitivity temperature drift	$1 \times 10^{-6/\circ}C \cdot Rnt$
Settling Time	1.5 s
Interfaces	RS-232C
Dimensions	263x487x322 mm (WxDxH)
Usable height of draft shield	235 mm
Size of weighing pan	78 x 73 mm

Standard Features:
SmartSens, hands-free operation SmartGrid, revolutionary grid weighing pan
SmartScreen, ease of use with color touch-screen display
ErgoClip for secure placement of all tare containers
LevelControl warns if balance is not levelled correctly Automatic and detachable draftshield proFACT
RS-232C interface
Option slot for second interface, choice of 7 different interfaces Inner draftshield plate, adjustable
Protective cover Prepared for weighing below the balance
Detachable and adjustable Terminal

All Models Offer:

QM-Toolbox, including user administration and password protection MinWeigh Target weighing Statistical evaluation Formulation with security check Density determination Percent weighing Piece counting with reference optimization Data entry of 4 ID's



IET | International Equipment Trading Ltd. <u>www.IetLtd.com</u> | *Proudly serving laboratories worldwide since 1979* CALL +847.913.0777 for Refurbished & Certified Lab Equipment

IET - Used Lab Equipment - Refurbished Analytical Laboratory Instruments